

Streamline all your stone management procedures.



Fusion[®] Quattro[®]

EXTRACTION BALLOON

DUAL PLATFORM ERCP

ACCESS

STONE MANAGEMENT

STRICTURE MANAGEMENT



Fusion® Quattro®

EXTRACTION BALLOON

Used for endoscopic removal of stones in the biliary system and for contrast injection.

- 4 inflation diameters in 1 device to meet your specific clinical needs
- Ease of injection with contrast exiting the distal tip for above injection balloons
- Radiopaque bands enhance fluoroscopic visualization of the extraction balloon
- Reference marks on syringe indicate each inflation diameter

Order Number	Reference Part Number	Balloon Diameter mm	Catheter Size Fr	Catheter Length cm	Injection Port	Wire Guide Diameter inch
G48225	FS-QEB-A	8.5/10/12/15	6.6	200	Above Balloon	.035
G48226	FS-QEB-B	8.5/10/12/15	6.6	200	Below Balloon	.035
G31921	FS-QEB-XL-A	12/15/18/20	6.6	200	Above Balloon	.035
G31922	FS-QEB-XL-B	12/15/18/20	6.6	200	Below Balloon	.035

Minimum accessory channel 3.2 mm. **Caution:** This product contains natural rubber latex, which may cause an allergic reaction.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com-A4



AI, ESC, IR, OHNS, PI, RH, SUR-A4